

sempermed®
derma PF

DISTRIBUTION PAR



STERIGENE
INGÉNIERIE DES PROCÉDÉS PROPRES & STÉRILES

DERMA PF FOR SENSITIVE HANDS

HIGHLY SKIN-FRIENDLY – the Semperit Post Production Process
EASY TO DON – patented synthetic inner coating
MAXIMUM COMFORT – fatigue-proof working



FOR SENSITIVE HANDS

The Sempermed® Derma PF offers good reasons to switch to powder-free gloves.

■ HIGHLY SKIN-FRIENDLY

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

■ EASY TO DON

With its synthetic inner coating the glove is easy to put on due to its patented net structure.

■ SURE GRIP

The Sempermed® Derma PF has more grip and thus combines safety, tactical properties and instrument guidance.

■ GOOD VALUE FOR MONEY

The long-term cost reduction is significant, if you take the very high follow-up costs of powder damage into consideration. The cost of purchasing Sempermed® Derma PF are minimal by comparison.



TECHNICAL DATA

sempermed® derma PF

TYPE	Sterile surgical glove for single use, powder-free, with synthetic lining	
MATERIAL	Natural rubber latex (NR)	
COLOUR	Natural white	
GLOVE SHAPE	Fully anatomical with rolled rim	
SIZE/OVERALL LENGTH as per EN 455-2	6 and 6½ 7, 7½ and 8 8½ and 9	min. 260 mm min. 270 mm min. 280 mm
WALL THICKNESS measured in single layer in the palm area	0.16–0.20 mm	
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 1.0 instead of AQL 1.5)	AQL 1.0	
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead ≥ 9 N)	≥ 12 N	
DURABILITY	3 years in original packaging if stored as per DIN 7716, ISO 2230	
STERILISATION	Radiation STERILE R	
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223	
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Directive 89/686/EEC	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾	
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical fibre-free peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs	
ARTICLE NUMBERS	SPTDERMAPF Tailles 6.0 / 6.5 / 7.0 / ... / 9.0	

CONTACT US!

 **STERIGENE**
INGÉNIERIE DES PROCÉDÉS PROPRES & STÉRILES

2 rue André Citroën
95130 FRANCONVILLE
Tél. : +33 (0)1 34 44 23 23
Fax : +33 (0)1 34 44 23 24
www.sterigene.com

¹⁾Please note the Sempermed Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock.**